

INCUWHITE

A More Natural Light for Incubators

Our engineers designed these white LEDs to match natural light inside our incubators. This creates an ideal environment for insects and other organisms. Compared to fluorescent lighting, the benefits of IncuWhite make our best-selling, reliable incubators an even more valuable investment.

Increased shelf space and growing height due to streamlined design

Significantly longer light lifetime


Brighter light intensity


More closely matches natural light

Flicker-free and dimmable in one-percent increments

More resistant to dust and water (protected from low-pressure jets)

Lower energy use



IncuWhite™ spectral output:**Installation****Incubators**

IncuWhite™ LEDs are the standard lighting system utilized in biological incubators.

Vertical Lighting

Vertical lighting (included on all incubator model numbers).

Horizontal Lighting

Horizontal lighting (included on all incubator model numbers) is affixed to the underside of the unit's wire shelving.

1. Unpackage each shelf and locate shelf hardware (pilaster clips)
2. Install pilaster clips:
 - a. Four clips, located at the same height are required to properly support each shelf